

# ULTRASONIC TENSO

## NEUTRAL BOOSTER DETERGENT

ULTRASONIC TENSO is specially formulated for use as an active detergent in cleaning tool by immersion in ultrasonic cleaning machines.

It is used as a booster washing product to increase the cavitation power of acid and alkaline products used in ultrasonic cleaning machines.

### DOSE AND WAY OF USE

ULTRASONIC TENSO must be used diluted in water ultrasonic tank at a proportion between 0.1 and 0.5% according to the tenacity of the soil to be removed.

### TECHNICAL SPECIFICATIONS

COLOUR: BROWN

APPEARANCE: TRANSPARENT LIQUID

SMELL: CHARACTERISTIC

DENSITY: 1,07 +/- 0,05 g/cm<sup>3</sup>

PH 1%: 7,2 +/- 0,5

### CONTAINS

15 % - 30 % Non ionic surfactants

Preservatives: Methyl isothiazolinone, Benzyl isothiazolone

This product meets the existing legislation on biodegradability and the Sanitary Technical Regulation on Detergents.

---

[www.tierratech.com](http://www.tierratech.com)